


true  
nxs

## SPINAL FIXATION SYSTEM

### Features & Benefits



Easiest system to use on the market with implants pre-loaded on instruments.

Superb biomechanical strength.

No contouring required.

Ultra-low profile.

Can add on to system in the event of adjacent level degeneration without removing system.

Simple, sterile, single use instrument kit providing precise tolerances for every case.

true™  
spinal fixation system



True NXS™ Rod



NXS™ Cap  
(on delivery tube)

START HERE



Basic Instrument Kit  
Min QTY: 1



Removal Instrument Kit  
Min QTY: 1



Pedicule Prep Probe  
Min QTY: 1

- Straight Lenke
- Curved Lenke
- Straight Paddle
- Curved Paddle
- MIS Pedicle Prep
- Jamshidi Needle



Cannulated Tap  
Min QTY: 1

- 5.5mm
- 6.5mm
- 7.5mm



Guide Wire  
Min QTY: 1

- Small
- Large
- XL



Rod Benders  
Min QTY: 1

true  
nxs

# CUSTOMER ORDER CHECKLIST

## NXS Trials



NXS Rods  
Min QTY: 2 Rods

true  
nxs



NXS Cap  
Min QTY:  
1 Cap / NXS Rod  
2 Caps / Starter Rod



NXS Starter Rods  
Min QTY: 2 Rods

## True Polyaxial Pedical Screw

### Available Sizes:



5.5mm x 25mm	6.5mm x 40mm
5.5mm x 30mm	6.5mm x 45mm
5.5mm x 35mm	6.5mm x 50mm
5.5mm x 40mm	7.5mm x 35mm
5.5mm x 45mm	7.5mm x 40mm
5.5mm x 50mm	7.5mm x 45mm
6.5mm x 35mm	7.5mm x 50mm



INNOVATIVE SURGICAL DESIGNS

(812) 369-4252 | info@innovativesd.net  
3903 S. Walnut St. | Bloomington, IN 47401